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
# Sakamoto I.V. Injection Trainer with Overswelling Vi-Pa

## User's Manual

Only use **water**\* or **artificial blood** for carrying out injection training using this product!  
Furthermore, do not use syringes with bent needle tips!

\* Use either tap water or distilled water for injection use.

### Table of contents

1	Introduction	22	6	How to use	28
2	Product contents	23	7	Cleaning and storage	37
3	Cautions on use	24	8	Troubleshooting	38
4	Product applications	26	9	Specifications	40
5	Part names	27			

SAKAMOTO MODEL  
CORPORATION

# 1 Introduction


## **Thank you very much for purchasing the " Sakamoto I.V. Injection Trainer with Overswelling Vi-Pa ".**

This user's manual contains cautions and notes on how to use this product safely. Please read this manual carefully in order to ensure that you know how to use the product safely.


In addition, after reading the manual, keep it in a safe place where it can be referred to at any time later on.

After purchasing the product, check the package contents carefully to make sure that nothing is missing or broken. If you notice any problems, contact Sakamoto Model.


## 2 Product contents


Arm model × 2


Overswelling injection pad × 2


Tank with integrated pump × 1

AC adapter × 1

AA-size batteries × 6

Vein and nerve sheet × 2


Cup × 1

Spare stopper cap × 3

### TIP

- When the vein and nerve sheet is placed over the arm model, the injection location can be checked. However, do not carry out injection while the sheet is still attached.

## 3 Cautions on use

In order to use this product safely, be sure to observe the following.

(Items which are particularly important are marked with "●".)

● **Do not subject the product to impacts or othe forces.**

- Do not drop.
- Do not bump against solid objects.
- Do not apply excessive loads.

If any of the above occur, damage may result.

● **Always use fresh needle tips in the syringes when practising.**

If you use a syringe with a needle that has a bent or damaged tip, it may cause the inside of the main unit to become damaged or the epidermis to become degraded.

● **Only use water or artificial blood as the fluid for injecting.**

If a fluid other than those mentioned above is used, damage or disagreeable odors may result. Use either tap water or distilled water for injection use. To create artificial blood, add 1 to 2 drops of the artificial blood concentrate to water.

● **Do not inject outside the injection area.**

The area where injection is possible is limited by the internal mechanism. If you inject anywhere outside the injection area, it may damage the product. (Refer to "Part names". (P.26))

● **Do not carry out injection while the vein and nerve sheet is still attached to the arm model.**

The injection needle cannot pierce the vein and nerve sheet. If the sheet is injected, it will damage the sheet and the syringe, and personal injury may also occur.

● **Do not disassemble or modify the product.**

If you attempt to disassemble or modify the product, it may become damaged or problems with operation may result. Never disassemble or modify the product.

● **Do not write on the arm model with items such as ballpoint pens.**

If you write on the arm model using a ballpoint pen or similar object, the surface of the arm model will become smeared with ink which cannot be removed.

● **Do not use continuously for long periods.**

As a means to protect the product, avoid using it for long periods of time. (No more than 3 hours at a time is recommended.)

- **Hold the body of the AC adapter when removing it from an outlet.**

When disconnecting the AC adapter from the wall outlet, hold it by the AC adapter body, not by the cable. If you hold by the cable when disconnecting it, the cable may become damaged, and it may lead to electric shocks or fire occurring. In addition, do not connect or disconnect the cable with wet hands.

- **If you notice any problems when using the product, stop using it.**

If problems occur with the product during use, stop using it, turn off the power, and disconnect the power cable from the electrical outlet. If the batteries have been used, remove them. Then contact Sakamoto Model. If you continue using the product, accidents may result.

- **Do not store in places which are exposed to high temperatures, high levels of humidity or large amounts of dust.**

If the product is stored in a place which is exposed to direct sunlight, near heating appliances or in places with large quantities of oil, moisture or dust, problems with operation may result.

- **There are limits to the effectiveness of simulator training.**

This product is a simulator capable of supporting advanced medical training. However, Sakamoto Model does not guarantee that training using this product will provide the same level of expertise as actual treatment.

- **Do not use this product for other applications.**

Do not use this product for applications other than the original application for which it is designed. The product may become damaged or injury may result.

## 4 Product applications

This product can be used to carry out the following training.


### ● Intravenous injections

### ● Drawing blood


- Overswelling can be created using the pump's negative pressure.
- The pump creates the effect of blood circulating through the veins, so that backflow can be checked when drawing blood.

## 5 Part names

### ■ Arm model


### ■ Overswelling injection pad

#### ⚠ CAUTION

Do not inject outside the injection area. Overswelling will not occur.

### ■ Tank with integrated pump

# 6 How to use

## ■ Preparation


### TIP

This product includes two arm models which can both be connected at the same time so that two people can carry out training simultaneously. The following describes how to connect a single arm model, but the connection procedure is the same for both.

### 1 Insert the batteries into the tank with integrated pump, or connect the AC adapter.


#### <When using the batteries>

Tank with integrated pump

#### <When using the AC adapter>


Tank with integrated pump


### ⚠ CAUTION

- Be sure to insert the batteries so that their polarities are correct.
- When using the batteries (alkaline) to operate two arms, the batteries will last for approximately 6 hours.
- Do not handle the AC adapter with wet hands. Otherwise electric shocks may occur.

- 2** Insert the connection tubes of the overswelling injection pad into the arm model.


- 3** Place the overswelling injection pad into the pad mounting hole and secure it with the fasteners.


**CAUTION**

Do not pull out the tubes on the underside of the overswelling injection pad. Otherwise it will damage the connections and air will leak from the pad, and overswelling will not be possible. However, if the tubes pull out, re-insert them as shown in the diagram.

- 4** Turn the stopper cap to remove it from the tank with integrated pump.


**CAUTION**

- Do not remove the stopper cap from the side which is not being used. Otherwise water will gush out when the pump is operated.
- Be careful not to lose the stopper cap.

**TIP**

- Three sockets are used for each arm model.
- When using only a single arm model at a time, either the left or right set of sockets can be used.
- By using two arm models, two people can carry out training simultaneously.

**5** Insert the connection tubes into the connection sockets.

**CAUTION**

- Insert the pink tube into the red socket, and insert the white tube into the black socket.
- Do not turn the connectors of the tubes any further than necessary.  
The tubes are made of plastic, and may become damaged.
- Always be sure to turn the tubes counterclockwise to unlock them first before disconnecting them. If you keep pulling them, the connectors will break.
- Connect the tubes carefully so that they do not break.


**6** Pour 200 ml of water into the cup.

**TIP**

- Use either tap water or distilled water for injection use. Artificial blood can also be used. To create artificial blood, add 1 to 2 drops of the artificial blood concentrate to water.
- In the following procedure, the term "water" is used exclusively.

**7** Pour the water into the water inlet of the tank with integrated pump, and then tighten the cap securely.


**CAUTION**

- Be careful not to spill the water. If the water gets spilled, wipe it away immediately. Otherwise electric shocks may occur.
- Tighten the cap securely. If the cap is not tightened securely, water may leak out.

## ■ Training


- 1 Turn on the power switch of the tank with integrated pump.

**TIP**


The pump will make a whirring noise while it is operating, but this is normal and not the sign of a malfunction.

- 2 Press the unit select button to select the number of arm models being used.

### <Overswelling the pad>


- 1** Press and hold the overswell button. ▶ The pad will gradually overswell.

**TIP**

- The amount of overswelling can be adjusted by how long the overswell button is held down. The pad will be fully overswollen after 20 to 30 seconds.
- This product is modeled on the left arm.


- 2** Press the overswell release button to release the overswelling of the pad.


## ■ Putting away

---


- 1 Turn off the power switch.


- 2 Set the tank with integrated pump so that the surface with the label is facing upward, and then turn on the power switch.


The water inside the connection tubes will be returned to the tank with integrated pump.

- 3** When all water has been removed from the tubes, turn off the power switch.

- 4** Open the cap of the tank with integrated pump and pour out all of the water.

**CAUTION**

Open the cap and remove all the water from the tank before pulling out the connection tubes. If the connection tubes are pulled out first, the water will gush out.

(A whooshing sound of air escaping will be heard when the cap is opened.)

- 5** From this point on, carry out steps **5** to **1** (in that order) in "Preparation".

## 7 Cleaning and storage

Be sure to observe the following when cleaning and storing the product.

- **Remove the batteries when not using the product.**

If the batteries leak inside the battery compartment, problems with operation will result.

- **If not using the product, disconnect the AC adapter from the wall outlet.**

When not using the product or when cleaning or storing it away, be sure to turn off the power switch and disconnect the AC adapter from the wall outlet.

- **Clean the arm models and the surface of the overswelling injection pad body by wiping them gently with a soft cloth.**

- If it is extremely dirty, wipe it with a cloth which has been moistened with diluted alcohol, and then wipe it again with a dry cloth.
- Do not use benzine, thinner or other solvents. Otherwise, it may damage the product or cause the epidermis to deteriorate.
- If you write on the product using a ballpoint pen or similar object, the surface of the product will become smeared with ink which cannot be removed.

- **Do not store in places which are exposed to high temperatures, high levels of humidity or large amounts of dust.**

If the product is stored in a place which is exposed to direct sunlight, near heating appliances or in places with large quantities of oil, moisture or dust, problems with operation may result.

- **The overswelling injection pad is a consumable part.**

The overswelling injection pad will deteriorate when it has been injected many times. Furthermore, the overswelling function will also decrease over time. In such cases, replace the overswelling injection pad with a new one (replacement part) which can be purchased separately.

## 8 Troubleshooting

No.	Symptom or problem	Remedy	Reference page
1	The product does not operate when the power switch is turned on.	Check that the batteries are inserted the correct way. (There may be no batteries inserted or the polarities may be incorrect.)	"6 How to use" (P.27)
		Check that the AC adapter is connected correctly to the tank with integrated pump and to the wall outlet.	"6 How to use" (P.27)
2	When the power switch is turned on, the indicator flashes yellow.	Replace the batteries.	"6 How to use" (P.27)
3	The pad does not overswell when the overswell button is pressed.	Check that the connection tubes are connected correctly and that they are securely installed.	"6 How to use" (P.30)
		Check if the tubes have pulled out from the underside of the overswelling injection pad.	"6 How to use" (P.28)
		Replace the overswelling injection pad with a new one (replacement part) which can be purchased separately.	"7 Cleaning and storage" (P.36)
4	Tubes have pulled out of overswelling injection pad.	Re-connect the connection tubes in the correct places. (Do not pull the pad connectors when reconnecting the tubes.)	"6 How to use" (P.28)

No.	Symptom or problem	Remedy	Reference page
5	Backflow does not occur.	Check that the connection tubes are not bent and that nothing is sitting on top of the tubes.	"6 How to use" (P.30)
		Check that the connection tubes are connected correctly and that they are securely installed.	"6 How to use" (P.30)
		If two arm models are connected, check that the unit select button is set to "2". (The indicator illuminates blue.)	"6 How to use" (P.32)
		Check that the cap of the tank with integrated pump is securely tightened.	"6 How to use" (P.31)
6	Water leaks out from the cap of the tank with integrated pump.	Check that the cap is securely tightened.	"6 How to use" (P.31)
7	There is a disagreeable odor.	Stop using the product and contact Sakamoto Model. (The inside of the main unit may be corroded.)	-

If any problems are not fixed when the remedies given above are carried out, contact Sakamoto Model.

## 9 Specifications

### ■ Overall product

Power supply	AC adapter	AC 100 – 240 V 50/60 Hz
	Batteries	DC 9 V
Power consumption		11 W
Weight		3.5 kg
Continuous operating time		3 hours or less recommended

### ■ Arm model

Dimensions	W 60 cm × D 11 cm × H 12 cm (Dimensions are approximate.)
Weight	840 g

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URL <http://sakamoto-model.co.jp>

E-Mail [trade@sakamoto-model.com](mailto:trade@sakamoto-model.com)

Product specifications and appearance are subject to change for improvement without notice.

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- Some of the components used in this product are manufactured by other companies, and product names, etc. for these components may be trademarks or registered trademarks of the manufacturers.