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The image shows a screenshot of a video call interface. At the top, there is a video feed of a woman with a microphone. Below the video feed, there is a control bar with three buttons: a mute button (a microphone with a slash), a cancel button (a telephone handset), and a speaker button (a speaker icon). The mute and speaker buttons are highlighted with dashed boxes. Two callout boxes point to these buttons. The left callout box says "カメラ側からの音声をOFFにします" (Turn off audio from the camera side). The right callout box says "カメラ側への音声をOFFにします" (Turn off audio to the camera side).

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"KOHARU Plus"

User's Manual

- Refer to the user's manual that is provided with the "M181 I.V.Pad" for further details on the accessory "M181 I.V.Pad".
- An SD card reader (sold separately) will be needed in order to transfer the recorded image data to a PC.

Table of contents

1	Introduction	41	6	How to use	52
2	Quick preparation guide.....	42	7	Cleaning and storage.....	65
3	Product contents	46	8	Troubleshooting	66
4	Cautions on use	48	9	Specifications	67
5	Product applications	51			

(Refer to the following page for a detailed Table of Contents.)

You can view this product's introductory page on our company's website by using either the following URL or the QR code* at right.

● Sakamoto Model Corporation's Introductory page for the "KOHARU Plus"
<https://sakamoto-model.com/product/nursing/m206/>

* QR code is a registered trademark of Denso Wave Inc.

Table of contents

1	Introduction	41
2	Quick preparation guide	42
3	Product contents	46
4	Cautions on use	48
5	Product applications	51
6	How to use	52
	6-1 Basic tablet operations	52
	6-2 Setting (switching) the camera images	53
	6-3 Adjusting the volume.....	56
	6-4 Playing back recorded data	58
7	Cleaning and storage	65
8	Troubleshooting	66
9	Specifications	67

1 Introduction

Thank you very much for purchasing the "KOHARU Plus".

This user's manual contains cautions and notes on how to use this product safely. Please read this manual carefully in order to ensure that you know how to use the product safely.

In addition, after reading the manual, keep it in a safe place where it can be referred to at any time later on.

After purchasing the product, check the package contents carefully to make sure that nothing is missing or broken. If you notice any problems, contact Sakamoto Model.

■ About this user's manual

● Cautions and other symbols

This manual uses the following symbols to indicate cautions on handling and other information which may be useful to know during use.

	This indicates a caution that applies when using the product. If the instructions listed are not followed, personal injury or damage to the product or other equipment may result.
---	--

	This indicates information which may be useful to know during use, or which should be taken note of in order to use the product properly.
---	---


● About illustrations

For reasons such as product improvements, the illustrations which appear in this manual may differ in some parts from the product which you have.

The details in the instructions will not have changed, so please follow the instructions.

2 Quick preparation guide


This guide describes the mechanism and how to start using the "KOHARU Plus".
Be sure to read the "Cautions on use" (P.48) before using this product.

TIP

- Using the objective camera is optional.
- If a howling noise occurs, connect the headset to the tablet.



2 Connect the objective camera power cable to the objective camera.

FAQ


- Wi-Fi does not connect →P.44
- Switching the camera →P.53
- Adjusting the volume →P.56
- Using recorded data for training →P.58

3 Connect the headset and turn on the power for the tablet.

Tap the "KO HARU Plus" application.

* Refer to the "KO HARU Plus" application login information card.

Voice Call
Tap "Voice Call" to start training!

< Wi-Fi connection is not enabled >


- 1 Check the wireless password which is written on the Wi-Fi information card that is included with the product.

TIP

The text written on the card may vary from the above diagram depending on changes to specifications.


- 2 Tap [Settings] on the tablet. ▶ Tap [Wi-Fi]. ▶ Tap the Wi-Fi switch to turn on Wi-Fi.

TIP

- For details on basic operation of the tablet, refer to "6-1 Basic tablet operations" (P.52).
- The screen displays may differ slightly depending on the specifications of the tablet and the OS.

3 Connect by following the procedure given below.


TIP

The "****" part varies for each product.


3 Product contents

If no quantity is listed, the quantity included is one.


Body

* The built-in main camera includes an SD card (32 GB microSD card).


Pajamas

Objective camera


* Includes an SD card (32 GB microSD card).


Tablet


Headset


Router power cable

Main camera power cable


Objective camera power cable


Tablet charging cable

SD card adapter ×2

M181 I.V. Pad


Silicone tape

Baby powder

Wi-Fi information card

"KOHARU Plus" application
login information card

TIP

- Refer to the instruction manual that is provided with the "M181 I.V.Pad" for the product contents for the "M181 I.V.Pad".
- In this manual, "SD card" is used as a general term which includes microSD cards.

4 Cautions on use

In order to use this product safely, be sure to observe the following.
(Items which are particularly important are marked with " ■ ".)

⚠ CAUTION

Be sure to read the cautions and notes which are given in the accessory instruction manual.


■ Body

The body has a skin which is made of silicone, and inside the body are several mechanical parts which include precision devices. Be sure to note the following during handling.


● Holding the body

In order to avoid damage, be sure to hold the body as shown below.

〈 When held by two people 〉

〈 When held by one person 〉

- Be careful not to apply excessive force to places such as the joints of the shoulders and thighs.
- Do not hold the body by the wrists or ankles when lifting it.

● Do not subject the body to severe impacts or other forces.

- Do not drop.
- Do not bump against solid objects.
- Do not apply excessive loads.

If any of the above occur, damage may result.

● Do not get the product wet or immerse it into water.

The body is not designed to be waterproof, so do not place it into water or get it wet, otherwise it will become damaged.

- **Do not write on the body with items such as ballpoint pens.**

If you write on the body using a ballpoint pen or similar object, the surface of the body will become smeared with ink which cannot be removed.

- **Do not touch against printed items.**

If the body is placed on top of or wrapped in printed items such as newspapers and books, the body may become imprinted with ink which cannot be removed.

■ Tablet

- **Do not use the tablet for other applications.**

If the tablet is used for other applications, problems may occur when using it to operate this product. In addition, do not install security software to the tablet or use it to connect to any Wi-Fi devices other than to this product's router.

- **Handle the tablet's screen carefully.**

Do not press the tablet's screen (touch panel) too firmly or use sharp objects such as fingernails or ballpoint pens to operate it. It will damage the tablet.

■ SD card

- **Do not apply too much force to the SD card.**

Do not place heavy objects on top of the SD card or bend it. Doing so might damage the model.

- **Avoid exposing the SD card to static electricity.**

If the SD card is exposed to static electricity, it may damage the SD card or the saved data may be erased.

- **The SD card are consumable items.**

Even if they are always used correctly, they may still become unusable due to normal wear and tear. If this happens purchase a suitable replacement, or contact Sakamoto Model.

■ General product handling

- **Do not process or modify the product.**

Do not make holes in the product, or process or modify it so that other items can be installed to it, etc. Damage or unexpected problems may occur.

- **Do not pull the cables with excessive force.**

When disconnecting the product's various cables, always be sure to hold the cables by the connectors, and pull them out gently. When disconnecting cables, do not hold the cables themselves, and do not apply excessive force to the cables, otherwise damage or problems with operation may occur.

- **Hold the body of the AC adapter when removing it from an outlet.**

When disconnecting the AC adapter from the wall outlet, hold it by the AC adapter body, not by the cable. If you hold by the cable when disconnecting it, the cable may become damaged, and it may lead to electric shocks or fire occurring. In addition, do not connect or disconnect the cable with wet hands.

- **Do not use the product if problems occur.**

If problems occur with the product during use, stop using it, turn off the power, and disconnect the power cable from the electrical outlet. Then contact Sakamoto Model. If you continue using the product, accidents may result.

- **Do not store in places which are exposed to high temperatures, high levels of humidity or large amounts of dust.**

If the product is stored in a place which is exposed to direct sunlight, near heating appliances or in places with large quantities of oil, moisture or dust, discoloration or deterioration may occur.

- **There are limits to the effectiveness of simulator training.**

This product is a simulator capable of supporting advanced medical training. However, Sakamoto Model does not guarantee that training using this product will provide the same level of expertise as actual treatment.

- **Do not use this product for other applications.**

Do not use this product for applications other than the original application for which it is designed. The product may become damaged or injury may result.

- **Be careful when handing recorded data.**

If there is a need for some reason to make recorded data public on the Internet, such as via a SNS, make sure that no invasions of privacy occur.

5 Product applications

The "KOHARU Plus" is a training model which allows communication to be carried out between the manikin and the tablet by means of audio and video, in order to practice communicating with patients.

● Simulated conversation training

By utilizing Wi-Fi communication, the instructor can pretend to be "KOHARU Plus" in order to carry out a conversation with the student.

● Image feedback training

Images which have been recorded using the camera which is built into the "KOHARU Plus" and the objective camera can be played back using the tablet or a PC in order to teach students about their manner of conversation.

● Practical nursing training

The following compatible products from Sakamoto Model can be attached to "KOHARU Plus" in order to carry out a variety of different types of practical training.

- M118-1 Intramuscular Injection Simulator "Musclemate" 2
- M148-3 I.V. Injection Pad "Veinmate" 2
- M153 Intramuscular Injection Simulator "Buttockmate"
- M153-1 Intramuscular Injection Simulator "Buttockmate" type 2
- M153-2 Buttocks Intramuscular Injection Simulator (wearing type) "Buttockmate" Right
- M156 Decubitus Treatment Simulator
- M169 Fit-on Manual Disimpaction Simulator
- M177 Postoperative Care Suit
- M178 Blood Pressure Simulator
- M181 Fit-on Intravenous Injection Trainer I.V.Pad
 - * Included with this product.
- M181-3 Fit-on Intravenous Injection Trainer I.V.Pad Mini
- M200-6 Female Catheterization & Enema Simulator Plus
- M200-7 Male Catheterization & Enema Simulator Plus

⚠ CAUTION

Do not attach any products or medical equipment other than the above compatible products.

● TIP

This product carries out remote conversations via Wi-Fi communication, so that a time lag of about 1 second occurs during audio conversations. Please bear this in mind when using the product.

6 How to use

6-1 Basic tablet operations

Operation of this product is carried out using the tablet. When you directly tap the buttons, etc. which are displayed on the tablet screen with your finger, the screen reacts so that you can enter the operation you would like to carry out.


This describes the basic terms and methods for operating the product using the tablet.

TIP

The display contents which are shown on the tablet screen are not displays which are used with this product. Use them for reference only.

(1) Tap


This is the most basic operation. Gently tap the tablet screen with your finger at the place where you would like to carry out an operation.

(2) Slide

Touch the tip of your finger against a control on the tablet screen, and then slide your finger while still pressing down.


This is used at times such as when adjusting the volume ("6-3 Adjusting the volume" (P.56) etc.).

(3) Pinch

Press two fingers against the tablet screen, and then pull them apart from each other (pinch out) or pull them in together (pinch in).

This is used at times such as when zooming in or out of the screen ("◆ Playing back using the "KOHARU Plus" application" (P.58)).

TIP


The camera images which are shown on the application screen in this chapter are sample images only. Actual camera images will be displayed during use.

6-2 Setting (switching) the camera images


At the time when this product is shipped from the factory, the images from the main camera and the objective camera are set to be displayed on the tablet. However, if they do not appear for some reason, carry out the following settings. This explains how to switch the camera settings.

The following explanation describes how to connect the objective camera. Carry out the same procedure for the main camera.



- 1 Tap the "  " (multi-screen) icon on the "KOHARU Plus" application screen.


- 2 Tap the "  " icon on the multi-screen.

3 Tap the camera that you wish to connect, and then select the check mark.


4 The tablet will connect to the camera and images will be displayed.

TIP

- If you tap "  " (single-screen display) at this time, only the selected camera images (the images surrounded by the blue lines) will be displayed.
- Conversation is carried out with the selected camera. Normally the main camera should be selected.
- If single-screen display has been selected, you can slide the side of the screen to switch cameras.

6-3 Adjusting the volume


◆ Adjusting the volume at the "KOHARU Plus" application screen

Tap "🔊" on the "KOHARU Plus" application screen to adjust the volume.

◆ Adjusting the volume during a conversation

Tap "🔊" on the "KOHARU Plus" application screen to adjust the volume.

TIP


- If you would like to turn off the sound coming from the camera while having a conversation with "KOHARU Plus", tap "🔇" (Mute). To turn off the sound going to the camera, tap "🔊" (Speaker). (Tap once more to cancel the previous operation.)

The sound coming from the camera will turn off

The sound going to the camera will turn off

- If the volume does not increase even after adjusting the volume of the microphone, adjust the volume of the tablet.

Use the side button to adjust the volume

- The headset dial can be used to adjust the volume.

Turn the dial to adjust the volume

6-4 Playing back recorded data

Playing back recorded data can be carried out using the "KOHARU Plus" application or a PC. Video playback software is required in order to play back the data using a PC.


TIP

The date and time cannot be set using the "KOHARU Plus" tablet. Therefore, note that the date and time display and the filename and updated date and time for the recorded files in the "KOHARU Plus" application will not be correct.

◆ Playing back using the "KOHARU Plus" application


The recorded data is recorded separately by both the main camera and the objective camera. Carry out the procedure in "6-2 Setting (switching) the camera images" (P.53) to switch the camera images.

- 1 Tap the "📁" (Playback & Memory) icon on the "KOHARU Plus" application screen.

2 Select the "Playback" tab. ▶ Search for the place to be played back in the playback screen.

(Playback will start automatically. Tap "⏸" (pause) to stop playback.)

TIP

- The date and time displayed on the screen will not be accurate.
- The recorded data will start being recorded automatically as soon as "KOHARU Plus" starts up. Once the SD card becomes full, the data will be overwritten in order starting from the oldest data.
Refer to "9 Specifications" (P.67) for the maximum recording time.

- 3 To return to the camera image screen, tap "<".

◆ Playing back using a PC

The recorded data which is saved in the main camera and the objective camera can be played back using a PC.

Follow the procedure given below to play back the recorded data. (The main camera is used as an example, but the procedure is the same for the objective camera.)

< Install "VLC". >

"VLC" is a free media player software which has been developed and distributed by VideoLAN which is specified as the software to use for playing back recorded data from "KOHARU Plus" on a PC. Follow the procedure given below to install the software.

- 1 Connect to the Internet, and download "VLC" from the following URL.
<https://www.videolan.org/vlc/>

⚠ CAUTION

The above URL is valid as of November 1, 2021.

- 2 Double-click the downloaded installer file to install the software. (Follow the instructions which appear on the screen to continue with the installation.)


TIP

The "*" portion of the filename will vary depending on the software version and the PC environment.

CAUTION


- Do not install "VLC" or any other video playback software applications onto the tablet.
- When "VLC" is started up for the first time, the following screen will be displayed. If you do not wish information to be collected by "VLC", deselect the "Allow metadata network access" check box. In the same way, deselect the "Regularly check for VLC updates" check box if you do not want updates.


< Playing back the recorded data on the PC >

The following explanation applies to the main camera. For the objective camera, carry out the procedure from step 3 onward.

1 Remove the chest cover from the body.


2 Put your hand inside the hole in the chest and pull out the main camera which is installed inside the neck.

TIP


There is no need to disconnect the cables which are connected to the main camera.

3 Take out the SD card from the camera.

4 Use an SD card reader or similar to read the SD card using the PC.


At PC screen

TIP

- If necessary, use the accessory SD card adapter.
- Use a separate commercially-available SD card reader.

5 Right-click on the file to be played back in order to select it. ▶ Click "Open with". ▶ Click "VLC".


6 The recorded file will then be played back.

TIP

- The recorded data is saved as a new file when "KOHARU Plus" is started up. If there are a large number of files remaining in the SD card, it can take time to find the file which you would like to view, so it is recommended that you regularly move or delete older files.
- When returning the SD card and the main camera to the body, carry out the steps in the order of ③ to ①. At this time, align the lens of the main camera with the hole in the neck of the body.

7 Cleaning and storage

CAUTION

Be sure to read the cautions and notes which are given in the accessory instruction manual.

■ Once training is finished

Be sure to observe the following points when leaving the product after training is finished. (Be sure to read the instructions in "When storing for long periods" also.)

● Disconnect the power cable.

When not using the product or when cleaning or storing it away, be sure to disconnect the power cable from the wall outlet.

● Turn off the tablet's power supply.

When not using the product or when cleaning or storing it away, turn off the power switch of the tablet.

■ When storing for long periods

When not using the product for long periods, clean it carefully and store it away while being sure to observe the following. (Be sure to read the instructions in "Once training is finished" also.)

● To clean the body, use a soft cloth to gently wipe the surface.

- If it is extremely dirty, wipe it with a cloth which has been moistened with diluted alcohol, and then wipe it again with a dry cloth.
- Do not use benzine, thinner or other solvents. In addition, do not apply alcohol directly to the surface. It may damage the product or cause the surface to deteriorate.
- If you write on the body using a ballpoint pen or similar object, the surface of the body will become smeared with ink which cannot be removed.

● Do not get the product wet or immerse it into water.

The body is not designed to be waterproof, so do not place it into water or get it wet, otherwise it will become damaged.

● Coat with baby powder to protect the body and to prevent stickiness.

If not using the product for a long period of time or if the body is sticky, coat it with the accessory baby powder. However, avoid applying baby powder to the camera lens hole in the neck of the body.

● Do not store in places which are exposed to high temperatures, high levels of humidity or large amounts of dust.

If the product is stored in a place which is exposed to direct sunlight, near heating appliances or in places with large quantities of oil, moisture or dust, discoloration or deterioration may occur.

● Lay the body sideways.

When storing the body, lay it sideways with all arm and leg joints fully extended. If the joints are bent, it may cause damage and deterioration.

8 Troubleshooting

No.	Symptom or problem	Remedy	Reference page
1	Camera images do not appear on the tablet. Sound cannot be heard from either the body or the tablet.	Are all of the cables (Camera power cable , Router power cable) connected correctly?	"2 Quick preparation guide"(P.42)
		Is the "KOHARU Plus" application running?	
		Check that the Wi-Fi connection for the tablet is enabled.	
		"Are the body or the objective camera too far away from the tablet, or is there an obstruction (such as metal or concrete) which is obstructing signal reception? Change the position of the body or objective camera and the tablet."	"9 Specifications"(P.67)
2	Main camera images do not appear on the tablet.	Is the position of the main camera lens aligned correctly with the hole in the neck of the body?	"Playing back the recorded data on the PC"(P.62)
3	No sound can be heard from either the body or the tablet.	Check the following audio adjustment statuses. <ul style="list-style-type: none"> • "KOHARU Plus" application • Tablet • Headset 	"6-3 Adjusting the volume" (P.56)
4	A howling noise occurs.	Connect the headset to the tablet.	"2 Quick preparation guide"(P.42)
5	SD card does not operate correctly. (The camera is not recognized, an error occurs, recording and playback are not possible.)	Try replacing the SD card. The SD cards are consumable items. When replacing an SD card, purchase a compatible commercially-available card and format it using the "KOHARU Plus" application. If you are not sure about the correct operation, please contact Sakamoto Model.	"9 Specifications"(P.67)

If any problems are not fixed when the remedies given above are carried out, contact Sakamoto Model.

9 Specifications

Body dimensions		Height approx. 160 cm
Body weight		17 kg
Recommended ambient temperature		0 – 40 °C
Power supply		AC 100 – 240 V 50/60 Hz
Power consumption	Cameras	5.4 W / 1 camera (10.8 W / 2 cameras)
	Wi-Fi router	10 W
	When charging tablet	20 W
Wi-Fi connection distance		Distance between tablet and router should be 5 m or less (this may vary due to obstacles or signal status).
SD card	Rating	The following microSD cards are compatible with this product. * microSDHC (8 GB to 32 GB) * microSDXC (64 GB to 128 GB)
	Format	FAT32
Maximum camera recording time	microSDHC (8 GB to 32 GB)	Maximum at high resolution approx. 96 hours * Approx. 3 hours per GB
	microSDXC (64 GB to 128 GB)	Maximum at high resolution approx. 384 hours * Approx. 3 hours per GB

TIP

Refer to the instruction manual that is provided with the "M181 I.V.Pad" for the specifications for the "M181 I.V.Pad".

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SAKAMOTO MODEL

C O R P O R A T I O N

Overseas Sales Department
1-9-29, Senbahigashi, Minoh City,
Osaka 562-0035, JAPAN

Tel +81-72-729-8350

Fax +81-72-729-8385

URL <https://www.sakamoto-model.com/>

E-Mail trade@sakamoto-model.com

Product specifications and appearance are subject to change for improvement without notice.

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- Some of the components used in this product are manufactured by other companies, and product names, etc. for these components may be trademarks or registered trademarks of the manufacturers.